

Jeunetique Cross-linked hyaluronic acid

Highly purified monophasic dermal fillers based on cross-linked hyaluronic acid





Contents

- 3** What is Meamo Labs Jeunetique Cross-linked hyaluronic acid?
- 5** The main advantages of Meamo Labs Jeunetique Cross-linked hyaluronic acid
- 9** Types of Meamo Labs Jeunetique Cross-linked hyaluronic acid
- 13** How to Store

What is Meamo Labs Jeunetique Cross-linked hyaluronic acid?




Meamo Labs Jeunetique cross-linked hyaluronic acid fillers are safe and effective, cutting-edge products designed to restore tissue volume.



Jeunetique dermal filler injections rehydrate soft tissues, restore their lost volume, reduce the appearance of wrinkles, and help create a defined, smooth facial contour.

Meamo Labs Jeunetique is a fully biocompatible and biodegradable dermal filler that can act over a long period, maintaining high uniformity.

The Meamo Labs Jeunetique line includes three products: Jeunetique Fine, Jeunetique Deep, and Jeunetique Ultra, which differ in the size of hyaluronic acid particles and their degree of cross-linking.



**The main advantages
of Meamo Labs
Jeunetique
Cross-linked
hyaluronic acid**



Optimal Viscoelastic Properties



Jeunetique is a monophasic dermal filler, which improves its homogeneity and viscoelasticity compared to similar biphasic products, offering several advantages:

- The filler is injected and distributed in the dermis as evenly as possible.
- Jeunetique is easily moldable, allowing for the achievement of natural-looking results.
- At the same time, the optimal viscoelasticity of the gel enhances its resistance to deformation under external forces.
- Additionally, the viscoelastic properties of the monophasic Jeunetique provide the maximum lifting effect.

Excellent Safety Profile



BDDE Scale

In the production of Jeunetique fillers, we have developed a new standard for removing unreacted residual materials, such as BDDE – the primary agent for cross-linking linear hyaluronic acid molecules.

We have extended the washing process at the final stage to 8 days, compared to the 1-2 days used by other manufacturers. As a result, the BDDE residue (<math><0.002\text{ mg per 1 ml}</math>) in the Jeunetique dermal filler line is significantly lower than the minimum permissible value, which leads to fewer side effects and a reduced risk of allergic reactions with our dermal fillers.




Long-lasting Effect

Jeunetique dermal filler line contains **24 mg of hyaluronic acid**, which promotes slower biodegradation compared to other fillers. The high concentration of the active ingredient ensures a longer-lasting result.



Painlessness and Comfort for Patients

The active formula of Meamo Labs Jeunetique fillers includes **0.3% lidocaine**, ensuring effective pain relief during injections, thus providing a more comfortable procedure.




Types of Meamo Labs Jeunetique Cross-linked hyaluronic acid



Meamo Labs presents
a professional line featuring
three products: Jeunetique
Fine, Deep, and Ultra

Jeunetique Fine

Jeunetique Fine is the filler with the lowest density in the Meamo Labs Jeunetique line, designed for the correction of delicate areas on the face and body. Its high plasticity allows for achieving natural-looking results.



CHARACTERISTICS

Concentration of Hyaluronic Acid:	24 mg/ml
Lidocaine:	0,3%
Packaging:	Single-use syringe (1.1 ml) , 2 × 30G needles
Injection Depth:	Intradermal, superficial-subdermal
Duration of Effect:	4-6 months

RECOMMENDED AREAS OF APPLICATION


- Lips
- Hands
- Marionette lines
- Nasolabial folds
- Tear troughs
- Medium-depth wrinkles

INDICATIONS

- Shallow wrinkles
- Crow's feet
- "Venus rings" on the neck
- Atrophic scars and stretch marks
- Skin laxity in the décolletage area
- Mild signs of aging
- Restoration of volume and hydration in the skin of the hands
- Marionette lines and expression lines around the mouth

Jeunetique Deep

Jeunetique Deep is a medium-density filler, ideal for intradermal and subdermal injection. It is highly effective for addressing more pronounced signs of aging, such as deep wrinkles and volume loss. Due to its plasticity, Jeunetique Deep perfectly meets the most common indications for correction, providing a long-lasting and natural result in areas requiring deeper treatment.



CHARACTERISTICS

Concentration of Hyaluronic Acid:	24 mg/ml
Lidocaine:	0,3%
Packaging:	Single-use syringe (1.1 ml), 2 × 27G needles
Injection Depth:	Intradermal, subdermal
Duration of Effect:	6-12 months

RECOMMENDED AREAS OF APPLICATION


- Nasolabial folds
- Marionette lines
- Temporal area
- Lips
- Mandibular line

INDICATIONS

- Physiological folds
- Lip volume enhancement
- Medium and deep wrinkles
- Moderate sagging of soft facial tissues

Jeunetique Ultra

Jeunetique Ultra is the highest-density filler in the Meamo Labs Jeunetique line, designed for medium and deep dermal implantation. With maximum cross-linking and larger particle size, Jeunetique Ultra delivers long-lasting results and superior stability, making it the ideal choice for addressing more pronounced signs of aging and for creating well-defined facial contours.



CHARACTERISTICS

Concentration of Hyaluronic Acid:	24 mg/ml
Lidocaine:	0,3%
Packaging:	Single-use syringe (1.1 ml) , 2 × 25G needles
Injection Depth:	Subdermal, supraperiosteal
Duration of Effect:	Up to 12 months

RECOMMENDED AREAS OF APPLICATION

- Middle third of the face
- Mandibular line
- Chin area
- Cheekbone area
- Angles of the lower jaw

INDICATIONS

- Pronounced sagging of soft facial tissues
- Deep age-related and constitutional changes
- Volume enhancement of the cheekbones, chin, and jawline

How to Store

Precautions for Safe Storage and Use of Meamo Labs Jeunetique Fillers

Keep away from heat sources and open flames, in a cool, well-ventilated area.



The filler should be used immediately after opening the sterile packaging.



The Meamo Labs Jeunetique injectable gel should be used before the expiration date indicated on the packaging.







Store at room temperature (up to 25°C / 77°F).

Do not freeze.





 www.meamo.space
 info@meamo.co
 [@meamolab](https://www.instagram.com/meamolab)
 +82 (70) 4513-1496